

D-72336 Balingen

Tel: +49-[0]7433-9933-0 Fax: +49-[0]7433-9933-149 E-Mail: info@kern-sohn.com Internet: www.kern-sohn.com

Operating instructions Baby scale

KERN MBC

MBC 15K2DNM MBC 20K10NM MBC 15K2DEM MBC 20K10EM

Version 4.1 2018-11 GB



MBC_NM-BA-e-1841

D	Weitere Sprachversionen finden Sie online unter www.kern-sohn.com/manuals
BG	Други езикови версии ще намерите в сайта <u>www.kern-sohn.com/manuals</u>
DK	Flere sprogudgaver findes på websiden <u>www.kern-sohn.com/manuals</u>
EST	Muud keeleversioonid leiate Te leheküljel www.kern-sohn.com/manuals
Е	Más versiones de idiomas se encuentran online bajo www.kern-sohn.com/manuals
GR	Άλλες γλωσσικές αποδόσεις θα βρείτε στην ιστοσελίδα <u>www.kern-sohn.com/manuals</u>
F	Vous trouverez d'autres versions de langue online sous www.kern-sohn.com/manuals
LV	Citas valodu versijas atradīsiet vietnē www.kern-sohn.com/manuals
FIN	Muut kieliversiot löytyvät osoitteesta www.kern-sohn.com/manuals
LT	Kitas kalbines versijas rasite svetainėje www.kern-sohn.com/manuals
GB	Further language versions you will find online under www.kern-sohn.com/manuals
RO	Alte versiuni lingvistice veți găți pe site-ul www.kern-sohn.com/manuals
I	Trovate altre versioni di lingue online in www.kern-sohn.com/manuals
SK	Iné jazykové verzie nájdete na stránke <u>www.kern-sohn.com/manuals</u>
NL	Bijkomende taalversies vindt u online op www.kern-sohn.com/manuals
SLO	Druge jezikovne različice na voljo na spletni strani <u>www.kern-sohn.com/manuals</u>
Ρ	Encontram-se online mais versões de línguas em <u>www.kern-sohn.com/manuals</u>
CZ	Jiné jazykové verze najdete na stránkách <u>www.kern-sohn.com/manuals</u>
PL	Inne wersje językowe znajdą Państwo na stronie <u>www.kern-sohn.com/manuals</u>
SE	Övriga språkversioner finns här: www.kern-sohn.com/manuals
н	A további nyelvi változatok a következő oldalon találhatók: www.kern-sohn.com/manuals
HR	Druge jezične verzije su dostupne na stranici: <u>www.kern-sohn.com/manuals</u>
NO	Andre språkversjoner finnes det på <u>www.kern-sohn.com/manuals</u>



KERN MBC

Version 4.1 2018-11 Operating instructions Baby scale

Contents

1	Technical data	. 5
2	Declaration of conformity	. 7
2.1	Explanation of the graphic symbols for medical devices	
3	Appliance overview	10
3.1	Overview of display	11
3.2	Keyboard overview	12
4	Basic Information (General)	13
4.1	Specific function	13
4.2	Proper use	13
4.3	Improper Use	14
4.4	Warranty	14
4.5	Monitoring of Test Resources	14
5	Basic Safety Precautions	15
5.1	Pay attention to the instructions in the Operation Manual	15
5.2	Personnel training	15
5.3	Preventing contamination	15
5.4	Appropriate use	15
6	Electromagnetic compatibility (EMC)	16
6.1	General hints	16
6.2	Electromagnetic interferences	17
6.3	Electromagnetic noise immunity	18
6.3.1	Crucial features of performance	20
6.4	Minimum distances	20
7	Transport and storage	21
7.1	Testing upon acceptance	21
7.2	Packaging / return transport	21
8	Unpacking, Setup and Commissioning	22
8.1	Installation Site, Location of Use	22
8.2	Unpacking	22
8.3	Scope of delivery	23
8.3.1	Modelle MBC-NM	
8.3.2	Modelle MBC-EM	-
8.4	Placing	
8.5	Rechargeable battery operation (MBC-A08)	
8.6	Battery operation	
8.7	Mains connection (MBC-NM)	27

8.8	Optional power supply units	. 27
8.9	Initial Commissioning	. 27
9	Operation	28
9.1	Weighing	. 28
9.2	Taring	. 28
9.3	Hold function (Standstill function)	. 29
9.4	Feeding function (control of weight gain)	. 29
9.5	Show another decimal place (not verified value)	. 30
9.6	Use the optional size measurement device MBC-A01	. 30
10	Menu	31
10.1	Navigation in the menu	. 31
10.2	Menu overview	. 32
11	Error messages	34
12	Service, maintenance, disposal	35
12.1	Cleaning	. 35
12.2	Cleaning/Disinfection	. 35
12.3	Sterilisation	. 35
12.4	Service, maintenance	. 35
12.5	Disposal	. 35
13	Instant help	36
14	Verification	37
14.1	Verification validity period (current status in G)	
15	Adjustment	39
16	Accessories (optional)	41

1 Technical data

KERN (Type)	MBC 15K2DNM	MBC 20K10NM	
Model	MBC 15K2DM	MBC 20K10M	
Weighing range (max)	6 kg / 15 kg	20 kg	
Minimum load (Min)	40 g / 100 g	200 g	
Readability (d)	2 g / 5 g	10 g	
Verification value (e)	2 g / 5 g	10 g	
Reproducibility	2 g / 5 g	10 g	
Linearity ±	2 g / 5 g	10 g	
Display	LCD with 25m	nm high digits	
Recommended adjustment weight, not added (class)	15 kg (M1)	20 kg (M1)	
Stabilization time (typical)	3 s	ec.	
Warm-up time	10 ו	min	
Operating temperature	10° C + 40° C		
Humidity of air	max. 80 % (not condensing)		
Input Voltage	100 V - 240 V, 50 / 60 Hz		
Auto Off	After "x" min adjustable without load change		
Bimensions fully mounted (W x D x H) mm890 x 470 x 175 (with integrated height measuring device)000 x 407 x 120 (without height measuring device)		nt measuring device) 07 x 120	
Dimensions display unit (W x D x H) mm	200 x 130 x 60		
Baby weighing pan (WxD) mm	600 x 280 x 55		
Weight kg (net)	4,6		
Rechargeable battery operation, optional	MBC-A08, internally 6x1.2 V 2000 mA		
Verification according to 2014/31/EU	Grade III		
Medical product in accordance with 93/42/EEC	Category I with measuring function		
Height measuring device, integrated	MBC-A01, Measuring range: 40 - 80 cm		

KERN (Type)	TMBC 15K2DEM-A	TMBC 20K10EM-A	
Model	MBC 15K2DEM	MBC 20K10EM	
Weighing range (max)	6 kg / 15 kg	20 kg	
Minimum load (Min)	40 g / 100 g	200 g	
Readability (d)	2 g / 5 g	10 g	
Verification value (e)	2 g / 5 g	10 g	
Reproducibility	2 g / 5 g	10 g	
Linearity ±	2 g / 5 g	10 g	
Display	LCD with 25m	nm high digits	
Recommended adjustment weight, not added (class)	15 kg (M1)	20 kg (M1)	
Stabilization time (typical)	3 s	ec.	
Warm-up time	10	min	
Operating temperature	10° C	. + 40° C	
Humidity of air	max. 80 % (no	ot condensing)	
Input Voltage	100 V - 240 V, 50 / 60 Hz		
Auto Off	After "x" min adjustable without load change		
Dimensions fully mounted (W x D x H) mm	890 x 470 x 175 (with integrated height measuring device) 600 x 407 x 120 (without height measuring device)		
Dimensions display unit (W x D x H) mm	200 x 130 x 60		
Baby weighing pan (WxD) mm	600 x 2	80 x 55	
Weight kg (net)	4,	6	
Rechargeable battery operation, optional	MBC-A08, internally 6x1.2 V 2000 mA		
Batteries	6 x 1.5 V AA		
Verification according to 2014/31/EU	Grade III		
Medical product in accordance with 93/42/EEC	Category I with measuring function		
Height measuring device, integrated	MBC-A01, Measuring range: 40 - 80 cm		

2 Declaration of conformity

To view the current EC/EU Declaration of Conformity go to:

www.kern-sohn.com/ce

The scope of delivery for calibrated weighing balances (= conformity-rated weighing balances) includes a Declaration of Conformity. Solely these weighing balances are classified as medical devices.

2.1 Explanation of the graphic symbols for medical devices



All medical scales marked in this way meet the requirements of the following directives:

- 1. 2014/31/EU: Directive on non-automatic weighing instruments
- 2. 93/42/EC: Directive concerning medical devices



Scales marked in this way underwent the conformity assessment procedure according to Directive 2014/31/EU for scales belonging to 3rd accuracy class.

WF 172795 Designation of the serial number of every device, applied at the device and on the packaging

(Number as an example)



Identification of the manufacturing date of the medical product.

(Year and month here as example)



"Please note the accompanying documents" or "Please note operating instructions"



"Please note operating instructions"



"Please note operating instructions"



Identification of manufacturer of medical product including address

Kern & Sohn GmbH D–72336 Balingen, Germany www.kern-sohn.com



"Electro-medical appliance" with attachment for type B



Device protection category II



Dispose of old appliances separately from your household waste!

Instead, take them to communal collection points.



Display of supply voltage for scales with polarity display.


Direct current supply voltage



Information




In order to avoid accidents, babies placed on the scales shall be continuously observed. Proceed in accordance with the indications specified on the scale!



Transportation protection




Power supply socket is marked with a small label on the side of the display




Level the scales before use

3 Appliance overview




- 1. Height measuring rod (optional)
- 2. Bubble level
- 3. Display Unit
- 4. LED
- 5. Baby weighing pan
- 6. Weighing pan
- Rubber feet (height adjustable)
 Mains connection (MBC-NM)
- 9. Battery compartment
- 10.RS232

3.1 Overview of display



Display	Description	Description
GROSS	Gross weight display	Lights up during indication of the gross weight of the baby (after drinking)
NET	Net weight display	Lights up during indication of the net weight of the baby (before drinking)
		Illuminated after weighing scale was tared
ZERO	Zeroing display	Should the balance not display exactly zero
		despite empty scale pan, press the →o↔ button. Your balance will be set to zero after a short standby time.
STABLE	Stability display	Scales are in a steady state
DRINK	DRINK function	Is displayed with active drink function
HOLD	HOLD function	Is displayed with active hold function
(D_)		Lights when the voltage drops below the prescribed minimum.
ű 🗆	Rechargeable battery symbol	Lights when the rechargeable battery capacity is almost exhausted.
		Lights when the rechargeable battery is fully charged.

3.2 Keyboard overview



MBC 15K2DNM MBC 20K10NM MBC 15K2DEM MBC 20K10EM

Key	Description	Function
	ON/OFF-switch	Turn on/off
→0←	Zeroing key	Weighing scale will be reset to "0.0" kg. For numeric entry: • Change decimal place
HOLD	HOLD button	Hold function
TARE	TARE button	Tare balance
[So)- CL	Feeding Function key	Differential weighing before and after the baby drinks
	20	The net weight of the baby will be shown: Before drinking
		In menu: Call up menu How to select menu items For numeric entry: Edit numeric value
		The gross weight is displayed: After drinking
		In menu:

⇒ Confirm selection

For numeric entry:

⇒ Confirm numerical value

4 Basic Information (General)



Weighing instruments have to be verified for the purposes stated below in accordance with Directive 2014/31/EU. Article 1, paragraph 4. "Determination of mass in the practice of medicine that is, weighing patients for reasons of medical supervision during medical surveillance, examination and treatment."

4.1 Specific function


- **Indication** Determining the body weight in the medical practice area.
 - Operated as "non-automatic weighing instrument" which means that you have to carefully put he baby in the centre of the weighing pan. Once a steady display value is shown, you can read the weight value.
- **Contra-** No contraindication known indication

4.2 Proper use

These scales serve as a means of determining the weight of babies in medical treatment rooms. The scales are suitable for recognising, preventing and controlling illnesses.



Scales fitted with a serial interface may only be connected to appliances in compliance with Directive EN60601-1.



4.3 Improper Use

Do not use these scales for dynamic weighing processes.

Do not leave permanent load on the weighing pan. This may damage the measuring system.

--

Impacts and overloading exceeding the stated maximum load (max) of the weighing plate, minus a possibly existing tare load, must be strictly avoided. This could cause damage to the balance.

Never operate balance in explosive environment. The serial version is not explosion protected. It should be noted that a flammable mixture of anaesthetics and oxygen or laughing gas may occur.

The structure of the balance may not be modified. This may lead to incorrect weighing results, safety-related faults and destruction of the balance.

The balance may only be used according to the described conditions. Other areas of use must be released by KERN in writing.

4.4 Warranty

Warranty claims shall be voided in case

- Our conditions in the operation manual are ignored
- The appliance is used outside the described uses
- The appliance is modified or opened
- mechanical damage and damage caused by media, liquids,
- natural wear and tear
- The appliance is improperly set up or incorrectly electrically connected
- The measuring system is overloaded
- Dropping the balance

4.5 Monitoring of Test Resources

In the framework of quality assurance the measuring-related weighing properties of the balance and, if applicable, the testing weight, must be checked regularly. The responsible user must define a suitable interval as well as type and scope of this test. Information is available on KERN's home page (<u>www.kern-sohn.com</u> with regard to the monitoring of balance test substances and the test weights required for this. In KERN's accredited DKD calibration laboratory test weights and balances may be calibrated (return to the national standard) fast and at moderate cost.

For balances with height measuring rods, we recommend a metrological examination of the accuracy of the height measuring rod, however, this is not mandatory as the determination of human body height involves rather large, intrinsic inaccuracies.

5 Basic Safety Precautions

⇔

5.1 Pay attention to the instructions in the Operation Manual



Carefully read this operation manual before setup and commissioning, even if you are already familiar with KERN balances.



5.2 Personnel training

The medical staff must apply and follow the operating instructions for proper use and care of the product.

5.3 Preventing contamination

To prevent cross-contamination (fungal skin infections, ...), clean the baby weighing pan or weighing platform every time after weighing. Recommendation: after a weighing procedure that could potentially result in contamination (e. g. after weighing that involves direct skin contact).

5.4 Appropriate use

- Before every use inspect the scales for any signs of damage.
- Maintenance and repeated verification Medical scales should be maintained and verified at regular intervals. (see Chapter 12.4)

6 Electromagnetic compatibility (EMC)

6.1 General hints



The installation and use of this baby scale MBC-NM requires special precautionary measures as outlined in the EMC information below.

This device complies with the limits set for medical electrical devices of group 1, class B (as per EN 60601-1-2).

Electromagnetic compatibility (EMC) describes a device's ability to perform reliably within an electromagnetic environment without causing inadmissible electromagnetic interference at the same time. Amongst other things, such disturbances may be emitted by connecting cables or the air.

Inadmissible disturbances from the environment may result in incorrect displays, inaccurate measured values or incorrect behaviour of the medical device. By the same token the medical device may in some cases cause such disturbances in other devices. To eliminate problems of that kind, we recommend you to take one or several of the measures listed below:

- Change the alignment or distance of the device to the source of EMI.
- Install or use the baby scale MBC-NM at a different location.
- Connect the baby scale MBC-NM to a different power source.
- For further questions please contact our customer services.

Disturbances may be caused by improper modification or add-ons to the device or not recommended accessories (such as power units or connecting cables). The manufacturer will not be responsible for these. Modifications may also result in a loss of authorisation relating to the use of the device.



Devices emitting high frequency signals (mobile telephones, radio transmitters, radio receivers) may cause interference in the medical device. For that reason do not use them near the medical device. Chapter 6.4 contains details about recommended minimum distances.

6.2 Electromagnetic interferences

Guidelines and manufacturer's declaration – electromagnetic interferences

The baby scale MBC-NM is designed for use in an electromagnetic environment that meets the requirements stated below. The customer or user of the medical electrical device must ensure that operation takes place in such an environment.

Emitted interference measurements	Conformity	Electromagnetic environment - guideline	
HF emissions as per CISPR 11 / EN 55011	Group 1	The baby scale MBC-NM uses HF energy merely for its internal working. Its HF emission therefore is very low and it is highly unlike to interfere with adjacent electronic devices.	
HF emissions as per CISPR 11 / EN 55011	Class B	The baby scale MBC-NM is designed for use in all equipment including those in living areas and those connected directly to the public power grid that also supplies buildings used for living purposes.	
Emission of harmonics as per IEC 61000-3-2	Class A		
Emission of voltage fluctuations / flicker	Conforms with		
as per IEC 61000-3-3			

Do not put the baby scale MBC-NM directly next to other devices and do not stack it with other devices. If this type of operation is necessary, observe the baby scale MBC-NM to ensure normal operation in such an arrangement.

6.3 Electromagnetic noise immunity

Guidelines and manufacturer's declaration - electromagnetic noise immunity

The baby scale MBC-NM is designed for use in an electromagnetic environment that meets the requirements stated below. The customer or user of the medical electrical device must ensure that operation takes place in such an environment.

Noise immunity tests	IEC 60601 test level	Conformity	Electromagnetic environment - guideline	
Discharge static electricity (DSE) as per IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV ± 8 kV	Floors should be made of wood or concrete or tiled with ceramic tiles. If floors are covered with synthetic material, relative air humidity must be at least 30%.	
Fast transient electrical disturbances / bursts as per IEC 61000-4-4	± 2 kV for power lines <u>+</u> 1 kV for input and output lines	± 2 kV <u>+</u> 1 kV	The quality of the supply voltage should match that of the typical business or hospital environment.	
Impulse voltages / surges as per IEC 61000-4-5	± 1 kV voltage Live wire - live wire ± 2 kV voltage Live wire - earth	± 1 kV Inapplicable	The quality of the supply voltage should match that of the typical business or hospital environment.	
Voltage dips, short-term disruptions and fluctuations in supply voltage as per IEC 61000-4-11	< 5 % U _T (> 95 % dip of U _T) for ½ period 40 % U _T (> 60 % dip of U _T) for 5 periods 70 % U _T (> 30 % dip of U _T) for 25 periods < 5 % U _T (> 95 % dip of U _T) for 5 s	Compliance with requirements under all postulated conditions Controlled switch off Return to undisturbed situation after user intervention.	The quality of the supply voltage should match that of the typical business or hospital environment. Where the user of the baby scale MBC-NM demands continuous operation even during disruptions to the power supply, we recommend powering the baby scale MBC- M by no-break power supply or battery.	
Magnetic field for supply frequency (50/60 Hz) as per IEC 61000-4-8	3 A/m	3 A/m 50/60 Hz	Magnetic fields for the supply frequency should match the typical values found in the particular business or hospital environment.	
NOTE U_T equals AC line voltage prior to application of test level.				

Guidelines and manufacturer's declaration - electromagnetic noise immunity

The baby scale MBC-NM is designed for use in an electromagnetic environment that meets the requirements stated below. The customer or user of the medical electrical device must ensure that operation takes place in such an environment.

Noise immunity tests	IEC 60601 test level	Conformity	Electromagnetic environment - guideline
Conducted HF disturbance variables as per IEC 61000-4-6	3 V _{mms} 150 kHz to 80 MHz	3 V	Do not use portable or mobile radio sets nearer to the baby scale MBC-NM or its wires than the distance recommended as safety distance which is calculated
Emitted HF disturbance variables as per IEC 61000-4-3	3 Vms 80 MHz to 2.5 GHz	3 V/m (((↔)))	according to the equation relevant for its transmission frequency. Recommended safety distance: $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ for 800 MHz to 2.5 GHz Use P as rated capacity of radio transmitter in Watt (W) as per details given by the radio transmitter manufacturer and d as recommended safety distance in metres (m). The field intensity of stationary radio transmitters should for all frequencies be lower according to an in situ ^a examination than the conformity level. ^b Interference may occur near devices bearing the symbol below.
 NOTE 1 Higher frequency range applies to 80 MHz and 800 MHz. NOTE 2 These guidelines may not be applicable in all cases. The spread of electromagnetic variables is influenced by absorption and reflections in buildings, objects and humans. ^a The field intensity of stationary radio transmitters such as base stations of wireless telephones and 			
 ^b For a frequency range of 150 kHz to 80 MHz field intensity should be less than 3 V/m. 			

6.3.1 Crucial features of performance



The baby scale MBC-NM does not have any crucial features of performance as per IEC 60601-1. The system may be subject to interference by other devices even if these devices conform to current emission requirements as per CISPR.

6.4 Minimum distances

Recommended safety distances between portable and mobile HF telecommunication devices and the medical device

The baby scale MBC-NM is designed for use in an electromagnetic environment in which HF disturbance variables are controlled. The customer or user of the medical electrical device can help avoiding electromagnetic disturbances by keeping the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the medical device – depending on the output performance of the communication device, as stated below.

Rated capacity of transmitter W	The safety distance depends on the transmission frequency m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00

For transmitters with a maximum rated capacity not stated in the table above you can calculate the recommended safety distance in metres (m) yourself by using the equation belonging to each column, whereby P equals the maximum rated capacity of the transmitter in Watt (W) as per details provided by the transmitter manufacturer.

NOTE 1 Higher frequency range applies to 80 MHz and 800 MHz.

NOTE 2 These guidelines may not be applicable in all cases.

The spread of electromagnetic variables is influenced by absorption and reflections in buildings, objects and humans.

7 Transport and storage

7.1 Testing upon acceptance

When receiving the appliance, please check packaging immediately, and the appliance itself when unpacking for possible visible damage.

7.2 Packaging / return transport



- ⇒ Keep all parts of the original packaging for a possibly required return.
- ⇒ Only use original packaging for returning.
- ⇒ Prior to dispatch disconnect all cables and remove loose/mobile parts.
- ⇒ Reattach possibly supplied transport securing devices.
- ⇒ Secure all parts such as the weighing pan, power unit etc. against shifting and damage.

8 Unpacking, Setup and Commissioning

8.1 Installation Site, Location of Use

The balances are designed in a way that reliable weighing results are achieved in common conditions of use. You will work accurately and fast, if you select the right location for your balance.

On the installation site observe the following:

- Place scales on a stable, even surface
- Avoid extreme heat as well as temperature fluctuation caused by installing next to a radiator or in the direct sunlight;
- Protect the balance against direct draughts due to open windows and doors;
- Avoid jarring during weighing;
- Protect the balance against high humidity, vapours and dust;
- Do not expose the device to extreme dampness for longer periods of time. Nonpermitted condensation (condensation of air humidity on the appliance) may occur if a cold appliance is taken to a considerably warmer environment. In this case, acclimatize the disconnected appliance for ca. 2 hours at room temperature.
- Avoid static charge of the balance and of the person to be weighed.
- Avoid contact with water.

Major display deviations (incorrect weighing results) may be experienced should electromagnetic fields (e.g. due to mobile phones or radio equipment), static electricity accumulations or instable power supply occur. In that case, the location must be changed.

8.2 Unpacking

Take the balance out of their packaging and place it at the intended position. When using the power pack, ensure that the power cable does not produce a risk of stumbling.

8.3 Scope of delivery

8.3.1 Modelle MBC-NM

- Balance
- Power pack unit (EN 60601-1 attestation of conformity)
- Operating instructions

8.3.2 Modelle MBC-EM

- Balance
- Batteries 6 x AA 1,5 V
- Operating instructions

8.4 Placing

Carefully remove the balance from the packaging, remove plastic cover and setup balance at the intended workstation.




Make sure that all transport locking devices are removed




To loosen the transport guard screw out transport screw [1] anticlockwise.

For transportation carefully screw-in transport screw clockwise till to the stopper and then fix it using locknut.



Levelling



Level balance with foot screws until the air bubble of the water balance is in the prescribed circle.

8.5 Rechargeable battery operation (MBC-A08)

(is possible by obtaining an optional battery power pack)



Open the battery compartment cover (1) at the base of the display unit and insert the rechargeable battery pack. Charge the battery for at least 12 hours before initial use.

The appearance of the symbol \square in the weight display indicates that the battery packs is almost exhausted. The weighing scale will remain ready for operation for a few more minutes before switching off in order to save battery. Load rechargeable battery.

Voltage has dropped below prescribed minimum.



Rechargeable battery very low.



Rechargeable battery completely reloaded

If the balance is not used for a longer time, take out the battery pack and store it separately. Leaking liquid could damage the balance.

8.6 Battery operation

As an alternative to rechargeable battery operation, the balance may also be operated with 6x AA batteries.

Open battery compartment cover (1) at the lower side of the display unit and insert batteries according to the example below. Lock the battery cover again. If the

batteries are empty, in the balance display appears the symbol batteries. To save battery power, the balance switches off automatically (see chap.11.6 Auto off).



Capacity of batteries exhausted.



Batteries will soon be flat.



Batteries are completely charged

Insert batteries

Remove battery compartment cover	
Connect battery holder to housing contact acc. to illustration	
Insert battery holder	
Insert batteries into battery compartment and lock with battery compartment	

cover.

8.7 Mains connection (MBC-NM)

Power is supplied by the external power unit which also serves to isolate the mains supply from the scale. The stated voltage value must be the same as the local voltage.

Only approved genuine KERN power supply units may be used in compliance with Directive EN 60601-1.

The small sticker attached to the side of the display unit indicates the power port:



The LED remains illuminated as long as the weighing scale remains connected to the mains.

The LED display provides information about the battery's charging status.

Green: battery is fully charged

Blue: battery is charging

8.8 Optional power supply units

Available power supply units (optionally)

- MBC-A04 (AUS/EU/UK/US/CH)
- MBC-A10 (EU/CH)

8.9 Initial Commissioning

In order to obtain exact results with the electronic balances, your balance must have reached the operating temperature (see warming up time chap. During this warming up time the balance must be connected to the power supply (mains, accumulator or battery) and be switched on.

The accuracy of the balance depends on the local acceleration of gravity. The value of gravity acceleration is shown on the type plate.


9 Operation

9.1 Weighing

STABLE ZERO GROSS 0.0000 kg	Start balance by pressing The balance will carry out a self-test The scales are ready for operation as soon as the weight display for "0.0kg" has appeared.
1	However, you can reset the weighing scale to zero by pressing the 40^{\pm} key.
	Put the baby in the centre of the weighing pan. Wait for stability display "STABLE", then read the weighing result.
1.	If the baby is heavier than the max. weighing range, the display shows "oL" (overload) and a beep sounds.

9.2 Taring

The tare weight of any preloads can be deducted by pressing a button so that the actual weight of the baby is displayed in subsequent weighings.




9.3 Hold function (Standstill function)

The balance has an integrated standstill function (mean value calculation). This allows one to weigh the baby exactly, even if it is not restful in the weighing pan.

9.4 Feeding function (control of weight gain)

The baby's weight can be saved before feeding. Then the weight gain can be calculated by pressing a button.







- \Rightarrow Place the baby on the weighing pan center before feeding.
- After the stability display shows STABLE", press [™]. The weight of the baby is recorded and stored. Display "DRINK" lights up.



- \Rightarrow Take the baby from the weighing pan.
- \Rightarrow Place the baby on the scale pan after feeding.
- ⇒ Press , the difference between the weight and the value before and after breastfeeding is displayed.




9.5 Show another decimal place (not verified value)

Press and hold for about 2 s whilst weighed result is being shown. The second decimal place will be shown for approx. 5 s.

9.6 Use the optional size measurement device MBC-A01

The scale has the ability to determine not only the weight but also the body height using the optional height measuring rod.



For this purpose proceed as follows:

- \Rightarrow Put the head stopper (left) (1) to zero (2)
- \Rightarrow Put the baby in the centre of the weighing pan.
- Move the height measuring rod (3) carefully to the right until the head stopper
 (3) gently touches the baby's head
- ⇒ With the right hand push the foot stopper (right) (4) carefully to the soles of the baby
- \Rightarrow On the scale read the baby's size.



For further information (for example, installation), refer to the instruction manual that comes with the height measurement.

10 Menu



Access to service menu "tCH"is locked in verified balances. To disable the access lock, destroy the seal and actuate the adjustment switch. For position of adjustment switch, see chap. 13.

Attention:

After destruction of the seal the weighing system must be reverified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.

10.1 Navigation in the menu		
Call up menu	⇒ Turn on the scale during the self-test press , the first function [F1 oFF] is displayed.	
Select function	⇒ With help of , the individual functions can be selected one after the other.	
Change settings	⇒ Confirm selected function by . The current setting wil be displayed.	
	⇒ Select the desired setting with $$ and press $$ to confirm or $$ to reject, the balance returns to the menu.	
Exit menu/ Return to weighing mode	\Rightarrow Press $\overline{\mathbf{T}}$, the balance will return to weighing mode.	

10.2 Menu overview

Function	Settings	Description	
	·	· ·	
F1 oFF	oFF 0*	Automatic shutdown off	
Automatic cutout	oFF 3	Automatic shutdown after 3 min	
Auto Off	oFF 5	Automatic shutdown after 5 min	
	oFF 15	Automatic shutdown after 15 min	
	oFF 30	Automatic shutdown after 30 min	
		·	
F2 bk	bl on	Back lighting for display on	
Background illumination	bl oFF	Display background illumination off	
of display	bl AU*	Backlighting for display will come on automatically as soon as the weighing scale is operated.	
tCH	D .	If display shows "Pin" adjust switch.	
Service menu	Pin	Then press , TARE, HOLD subsequently.	
P1 Spd	15*		
Display speed	30	Not documented	
	60		
	7.5		

P2 CAL	duA in	dESC	C 0.00
			C 0.000
			C 0.0000
			C 0
			C 0.0
		inC	Sd iv 1 div 1, 2, 5, 10, 20, 50
			Sd iv 2 div 1, 2, 5, 10, 20, 50
		CAP	CAP 1
		0/ 1	CAP 2
		CAL	UnLoAd
		StrAnG	St 100
		OUAIO	St 200
			St 500
	duA rA	dESC	C 0.00
	du/(I/(uroo	C 0.000
			C 0.0000
			C 0
			C 0.0
		inC	Sd iv 1 div 1, 2, 5, 10, 20, 50
			Sd iv 2 div 1, 2, 5, 10, 20, 50
		CAP	CAP 1
		0/ 1	CAP 2
		CAL	UnLoAd
		Justierung	
		StrAnG	St 100
			St 200
			St 500
	SnG rA	dESC	C 0.00
			C 0.000
			C 0.0000
			C 0
			C 0.0
		inC	Sd iv 1 div 1, 2, 5, 10, 20, 50
			Sd iv 2 div 1, 2, 5, 10, 20, 50
		CAP	CAP 1
			CAP 2
		CAL	UnLoAd
		StrAnG	St 100
			St 200
			St 500

P3 Pro	tri	Not documented
	CoUnt	Not documented
	rESEt	Reset weighing scale to factory setting
	SEtGrA	Not documented

* default setting

11 Error messages

Display	Description
Erry	 Zero range exceeded (on start-up or when pressing the key) Load on weighing pan Excess load, during zero setting of weighing scale Incorrect adjusting process Fault on load cell
Err8	 Value outside the A/D changer range Damaged weighing cell Damaged electronics
Err 19	 Unable to initialise zero point Measuring cell defective / overloaded Object on weighing pan / contact Transport safety device has not been removed Main board defective

Should other error messages occur, switch balance off and then on again. If the error message remains inform manufacturer.

12 Service, maintenance, disposal

12.1 Cleaning



Before any maintenance, cleaning and repair work disconnect the appliance from the operating voltage.

12.2 Cleaning/Disinfection

The scales plate (e.g. stool) and the housing shall be cleaned only by means of a cleaning agent for household use or a disinfecting agent available on the market, e.g. 70% solution of isopropanol. We recommend using the disinfecting agent intended for disinfection by means of wet wiping the surface. Observe the manufacturer's recommendations.

Do not use any polishing or aggressive cleaning agents such as vodka, white spirit or similar as they may damage the high quality surface.

In order to avoid cross contamination (mycosis), it is necessary to observe the following disinfection schedule:

- the scales plate before and after each measurement with direct contact with skin
- As necessary:
 - o display
 - o foil-covered keyboard



Do not spray any disinfecting agents on the device.

The disinfecting agent cannot get into the scales interior. Immediately remove contamination.

12.3 Sterilisation

It is not allowed to sterilize the device.

12.4 Service, maintenance

The appliance may only be opened by trained service technicians who are authorized by KERN.

We recommend regular verification of conformity with technical safety requirements (STK).

Disconnect the scales before opening.

12.5 Disposal

Disposal of packaging and appliance must be carried out by operator according to valid national or regional law of the location where the appliance is used.

13 Instant help

In case of a fault in the program sequence, the balance should be shortly switched off. The weighing process must then be restarted from the beginning.

Failure:	Possible causes:
The displayed weight does not glow.	 The balance is not switched on. The mains supply connection has been interrupted (mains cable not plugged in/faulty). Power supply interrupted. Rechargeable battery/batteries inserted incorrectly or empty No rechargeable battery/batteries inserted
The displayed weight is permanently changing	 Draught/air movement Table/floor vibrations The weighing pan is in contact with foreign bodies or is not correctly positioned. Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)
The weighing result is obviously incorrect	 The display of the balance is not at zero. Adjustment is no longer correct. Great fluctuations in temperature. The balance is on an uneven surface. Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)

Should other error messages occur, switch balance off and then on again. If the error message remains inform manufacturer.

14 Verification

General introduction:

According to EU directive 2014/31/EU balances must be officially verified if they are used as follows (legally controlled area):

- a) For commercial transactions if the price of goods is determined by weighing.
- b) For the production of medicines in pharmacies as well as for analyses in the medical and pharmaceutical laboratory.
- c) For official purposes
- d) For manufacturing final packages
- e) Determination of mass in the practice of medicine that is, weighing patients for reasons of medical supervision during medical surveillance, examination and treatment.

In cases of doubt, please contact your local trade in standard.

Verification notes:

An EU type approval exists for balances described in their technical data as verifiable. If a balance is used where obligation to verify exists as described above, it must be verified and re-verified at regular intervals.

Re-verification of a balance is carried out according to the respective national regulations. For verification validity period, s. chap. 16.1.

The legal regulation of the country where the balance is used must be observed!




Verification of the balance is invalid without the seal.

The seal marks attached on balances with type approval point out that the balance may only be opened and serviced by trained and authorised specialist staff. If the seal mark is destroyed, verification loses its validity. Please observe all national laws and legal regulations. In Germany a re-verification will be necessary.

Balances with obligation to verify must be taken out of operation if:

- The weighing result of the balance is outside the error limit. Therefore, in regular intervals load balance with known test weight (ca. 1/3 of the max. load) and compare with displayed value.
- The reverification deadline has been exceeded.

Position adjustment switch and seals:



- 1. Self-destroying seal mark
- 2. Cover
- 3. Adjustment switch

14.1 Verification validity period (current status in G)

Personal scales (including chair and wheelchair scales) in hospitals	4 year
Personal scales, when not located in hospitals (for example, doctor's offices and nursing homes)	unlimited
Baby weighing scales and mechanical birth weight scales	4 year
Bed scales	2 year
Scales in dialysis stations	unlimited

Rehab clinics and health authorities are treated as hospitals. (4 years of verification validity)

Not treated as hospitals (verification validity not limited) are dialysis stations, nursing homes and doctor's surgeries.

(Details derived from: "Information by the verification authority, weighing scales applied in medical use")

15 Adjustment

As the acceleration value due to gravity is not the same at every location on earth, each display unit with connected weighing plate must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the weighing system has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the display unit periodically in weighing operation.


1	•	Prepare the required adjustment weight. The adjustment weight to be applied depends on the capacity of a weighing scale, see chap. 1. Carry out adjustment as closely as possible to admissible maximum load of weighing scale. Information about test weights you will find in the internet under <u>http://www.kern-sohn.com</u>
	•	Observe stable environmental conditions. For warm-up time required for stabilisation see chpt 1.



Access to service menu "tCH"is locked in verified balances. To disable the access lock, destroy the seal and actuate the adjustment switch. Position of the adjustment switch see chap. 13. **Attention:**

After destruction of the seal the weighing system must be re-verified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.

Procedure:



	⇒ Ensure that there are no objects on the weighing pan.
UnloAd	⇔ Wait for stability display "STABLE", then confirm with
Ŭ 0.00.00.₀	⇒ The size of the currently set adjustment weight is displayed, the active site flashes.
(example)	If required, select with the digit to be altered and change the digit with .
	Confirm with , [LoAd] appears
LoRd	Put the required adjustment weight carefully in the centre of the weighing pan.
•	⇒ Wait until stability display "STABLE" appears
PRSS	⇔ Confirm with [PASS] is displayed.
GROSS	After the adjustment the balance will carry out a self-test. Remove adjusting weight during selftest, balance will return into weighing mode automatically.
	An adjusting error or incorrect adjusting weight will be indicated by the error message; repeat adjustment procedure.
	An adjusting error or incorrect adjustment weight will generate an error message ("Err 4"), repeat the adjustment process.

16 Accessories (optional)

Item number	Product
MBC-A01	Height rod
MBC-A04	Power supply unit (EU/UK/CH/US/AUS)
MBC-A05	Foot/head rest
MBC-A08	Rechargeable battery
MBC-A10	Power supply unit (EU/CH)