



USER MANUAL MS3910 Stand-on Floor Scale

 \fbox Please keep the instruction manual at hand all the time for future reference.

Explanation of Graphic Symbols on Label/Packaging

\triangle	Caution, consult accompanying documents before use	X	Separate collection for waste of electrical and electronic equipment, in accordance with Directive 2002/96/EC
	Manufacturer of medical device		Manufacturing year of medical device
	Carefully read user manual before installation and usage, and follow instructions for use.	木	Medical electrical equipment with Type B applied part
REF	Device catalogue number	EC REP	Authorized representative in the European Community
LOT	Manufacturer's batch or lot number	MD	Device is a medical device
SN	Serial number	UDI	Unique Device Identifier
	CE 2460		93/42/EEC as amended ical Device Directive. Four to Notified Body.
		Device complies with International Organization of Legal Metrology (Class III) requirements (verified models only)	
	/180122	Device complies with models only)	h EC directives (verified
		M: Conformity label Directive 2014/31/E weighing instrumen	U for non-automatic
			onformity verification was CE label was applied. (ex:
		0122: Refers to Not	tified Body for metrology

Copyright Notice Charder Electronic Co., Ltd.

No.103, Guozhong Rd., Dali Dist., Taichung City 41262 Taiwan Tel: +886-4-2406 3766 Fax: +886-4-2406 5612 Website: www.chardermedical.com E-mail: info_cec@charder.com.tw

Copyright[®] Charder Electronic Co., Ltd. All rights reserved. This user manual is protected by international copyright law. All content is licensed, and usage is subject to written authorization from Charder Electronic Co., Ltd. (hereinafter Charder) Charder is not liable for any damage caused by a failure to adhere to requirements stated in this manual. Charder reserves the right to correct misprints in the manual without prior notice, and modify the exterior of the device for quality purposes without customer consent.



Charder Electronic Co., Ltd. No. 103, Guozhong Rd., Dali Dist., Taichung City, 41262 Taiwan

CONTENTS

I. Safety Notes	5
A. General Information	
B. EMC Guidance and Manufacturer's Declaration	
II. Installation	
A. Assembly	
B. Inserting Batteries	
C. Using Adapter	
D. Attaching Height Rod to Column E. Attaching Thermal Printer	
III. Indicator	
A. Indicator and Key Functions	
B. Display layout	
IV. Using Device	
A. Basic Operation	
B. Hold	.23
С. ВМІ	
D. Tare	
E. Pre-Tare	
F. Print	
V. Device Setup	
B. Device Setup	
VI. Setup USB Connection to PC	
VII. Wireless Connection	
VIII. Troubleshooting	
IX. Product Specifications	
X. Declaration of Conformity	.40



A. General Information

Thank you for choosing this Charder Medical device. It is designed to be easy and straightforward to operate, but if you encounter any problems not addressed in this manual, please contact your local Charder service partner.

Before beginning operation of the device, please read this user manual carefully, and keep it in a safe place for reference. It contains important instructions regarding installation, proper usage, and maintenance.

Intended Use

This device is intended to measure the weight of subjects who can stand unassisted, for diagnosis of weight-related issues by professionals.

General Handling

- Device should be placed on stable, flat, solid, non-slippery surface.
- Usage on soft surfaces (ex: carpet) may result in inaccurate results.
- Ensure all parts are properly locked and tightened before operating the device.
- Device is intended to measure one subject at a time.

Safety Instructions

- Batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Expected service life: 5 years.
- Always comply with appropriate regulations when using electrical components under increased safety requirements.
- Ensure voltage marked on power supply matches mains power supply.
- The device is intended for indoor use only.
- Observe permissible ambient temperatures for use

Environmental

 All batteries contain toxic compounds; batteries should be disposed of via designated competent organizations. Batteries should not be incinerated.

Cleaning

Device surface should be cleaned using alcohol-based wipes. Corrosive cleansing liquids should not be used. Pressure-washers should not be used.

- Do not use large amounts of water when cleaning the device, as it may cause damage to the internal electronics.
- Always disconnect device from mains power before cleaning.

Maintenance

Device does not require routine maintenance. However, regular checking of accuracy is recommended; frequency to be determined by level of use and state of device. If results are inaccurate, please contact local distributor.

Warranty/Liability

- The period of warranty shall be eighteen (18) months, beginning on the date of purchase. Please retain your receipt as proof of purchase.
- No responsibility shall be accepted for damage caused through any of the following reasons: unsuitable or improper storage or use, incorrect installation or commissioning by the owner or third parties, natural wear and tear, changes or modifications, incorrect or negligent handling, chemical, electrochemical, or electrical interference.
- All maintenance, technical inspections, and repairs should be conducted by an authorized Charder service partner, using original Charder accessories and spare parts. Charder is not liable for any damages arising from improper maintenance or usage.

Disposal

This product is not to be treated as regular household waste, but should be taken to a designated collection points for electronics. Further information should be provided by local waste disposal authorities.

≜ Marning

- Only the original adapter should be used with the device. Using an adapter other than the one provided by Charder may cause malfunction.
- Do not touch the power supply with wet hands.
- Do not crimp the power cable, and avoid sharp edges.
- Do not overload extension cables connected to the device.
- Route cables carefully, to avoid tripping.
- Keep device away from liquids.
- Do not remove the plug by yanking on the cable.
- Use only a correctly wired (100-240VAC) outlet, and do not use a multiple outlet extension cable.

- Do not under any circumstances dismantle or alter the device, as this could result in electric shock or injury as well as adversely affect the precision of measurements.
- Do not place the device in direct sunlight, or in close proximity to an intense heat source. Excessively high temperatures may damage the internal electronics.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, EU representative (if device is used in EU member state), and competent authority of user/subject's member state.

B. EMC Guidance and Manufacturer's Declaration

Guidance and manufacturer's declaration-electromagnetic emissions

The MS3910 Stand-on Floor Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance	
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly	
Harmonic emissions IEC 61000-3-2	Class A	connected to the public low-voltage power supply network that supplies buildings used for domestic	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Compliance	buildings used for domestic purposes.	

Guidance and manufacturer's declaration-electromagnetic immunity

The MS3910 Stand-on Floor Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601	Compliance	Electromagnetic
	test level	level	environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	<u>± 8 kV contact</u> <u>± 2 kV, ± 4 kV,</u> <u>± 8 kV, ± 15 kV</u> <u>air</u>	<u>± 8 kV contact</u> <u>± 2 kV, ± 4 kV,</u> <u>± 8 kV, ± 15 kV</u> <u>air</u>	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%

Electrical fast transient/burst IEC 61000-4-4± 2kV for power supply lines+ 2kV for power supply lines + 1kV for input/output linesMains power quality should be that of a typical commercial or hospital environment.Surge IEC 61000-4-5± 1kV line(s) to line(s)+ 1kV line(s) to line(s)Mains power quality should be that of a typical commercial or hospital environment.Surge IEC 61000-4-5± 1kV line(s) to line(s)+ 1kV line(s) to line(s)Mains power quality should be that of a typical commercial or hospital environment.Voltage Dips, short interruptions and voltage power supply0% UT for 0.5 0% UT for 1 cycle0% UT for 1 0% UT for 1 cycle0% UT for 1 cycleMains power quality should be that of a typical commercial or hospital environment.						
	transient/burst	power supply lines + 1kV for input/output	supply lines + 1kV for input/output	that of a typical commercial or		
short interruptionscyclecyclethat of a typical commercial orand voltage0% UT for 10% UT for 1hospital environment. If thevariations oncyclecycleuser of the device requires		line(s) ± 2kV line(s) to	line(s) + 2kV line(s) to	that of a typical commercial or		
input lines IEC70% UT(30%70% UT(30%power mains interruptions, it is61000-4-11dip in UT) for 25 cyclesdip in UT) for 25 cyclesrecommended that the device be powered from an uninterruptible power supply or a battery.0% UT for 5 s0% UT for 5 sa battery.	short interruptions and voltage variations on power supply input lines IEC	cycle 0% UT for 1 cycle 70% UT(30% dip in UT) for 25 cycles	cycle 0% UT for 1 cycle 70% UT(30% dip in UT) for 25 cycles	that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or		
Power frequency(50/60 Hz) magnetic field30 A/mThe device power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.NOTE UT is the a.c. mains voltage prior to application of the test level.	frequency(50/60 Hz) magnetic field IEC 61000-4-8			magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.		

IOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration-electromagnetic immunity

The MS3910 Stand-on Floor Scale is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that is used in such an environment.

Immunity test	IEC 60601 test	Compliance	Electromagnetic
	level level		environment-guidance
Conducted RF	3 Vrms	3 Vrms	Portable and mobile RF
IEC 61000-4-6	150 KHz to 80 MHz	150 KHz to 80 MHz	communications equipment
Radiated RF IEC	<u>6 V in ISM bands</u>	IVITIZ	should be used no closer to any
61000-4-3	between 0,15 MHz	6 V in ISM	part of the device including
	and 80 MHz	bands between	cables, than the recommended
	<u>80 % AM at 1 kHz</u>	<u>0,15 MHz and</u>	separation distance calculated
	3 V/m	<u>80 MHz_</u> 80 % AM at 1_	from the equation applicable to

	80MHz to 2,7 GHz	<u>kHz</u>	the frequency of the
		<u></u>	transmitter.
		3 V/m	
		<u>80MHz to 2,7</u> GHz	Recommended separation
			distance:
			$d = 1,2 \sqrt{P}$
			$d = 1,2 \sqrt{P}$ 80MHz to 800 MHz
			$d = 2,3 \sqrt{P}$ 800MHz to 2,5 GHz
			Where P is the maximum output
			power rating of the transmitter
			in watts (W) according to the
			transmitter manufacturer and ${\it d}$
			is the recommended separation
			distance in metres (m).
			Field strongths from fixed DE
			Field strengths from fixed RF transmitters, as determined by
			an electromagnetic site survey ^a ,
			should be less than the
			compliance level in each
			frequency range ^b .
			Interference may occur in the
			vicinity of equipment marked
			with the following symbol:
			((t,y))
	Iz and 800 MHz, the h		
0	ption and reflection from	2	Electromagnetic propagation is
-	•		ations for radio (cellular/cordless)
-			and FM radio broadcast and TV
			cy. To assess the electromagnetic

considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and the MS3910 Stand-on Floor Scale

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of	Separation distance according to frequency of transmitter m				
transmitter W	150 kHz to 80 MHz d =1,2√P	80 MHz to 800 MHz d =1,2√P	800 MHz to 2,5 GHz d =2,3√P		
0,01	0,12	0,12	0,23		
0,1	0,38	0,38	0,73		
1	1,2	1,2	2,3		
10	3,8	3,8	7,3		
100	12	12	23		

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


II. Installation

A. Assembly

1. Attach column to base.




2. Secure column and base using screws




B. Inserting Batteries


1. Open battery housing cover




2. Accessing batteries




3. Use either rechargeable battery pack, or AA batteries




4. Ensure batteries are installed into the housing correctly



5. Install the battery housing into the compartment, and make sure the right side of housing pin is facing towards inside of the connecting position



6. Slide back the cover to close the battery housing compartment. Turn on power to confirm that battery is correctly installed.




Using Rechargeable Battery (optional)

The rechargeable battery should be recharged at least once every 3 months, regardless of if the device has been used. Battery can be charged by plugging device's exclusive adapter into AC Connector Port.

After a long period in storage (e.g. >3 months), the battery should run a full cycle (charge/discharge) to allow it to restore full capacity.

Ensure rechargeable battery housing is installed and inserted properly into the compartment.




If **Lo** prompt displays on the LCD, please charge battery promptly to avoid battery damage.


C. Using Adapter


Connect adapter to indicator before connecting to mains power supply
Disconnect adapter from mains power supply before unplugging

adapter pin from indicator.



D. Attaching Height Rod to Column






1. Attach brackets to column with round-head screws


2. Attach height rod to brackets using flat-head screws

Item	Name	Quantity
1	M5x0.8x11 round head screw	4
2	Relief Bushing	2
3	Bracket for HM200D/HM201D/HM201M	2
4	Height Rod (Compatible with: HM200D/HM201D/HM201M)	1
5	M5x10L flat head screw	2
6	Fixing block	2
7	M5x0.8x11	4

Connecting digital height rod to indicator (HM200D/HM201D)




E. Attaching Thermal Printer



Item	Parts	Qty
1	M5*15L head screw	1
2	Screws for printer bracket	2
3	Printer bracket	1
4	TP2100/TP2110 Thermal Printer	1 (purchased separately)



2. Install the thermal printer on the bracket



III. Indicator


A. Indicator and Key Functions

POWER PRE- TARE HOLD NET WEIGHT	chard	er	
		BMI	HOLD
	ZERO	2	3
חרו	M I-5 4	5	6
170. 5 ► cm	PRE- TARE 7	8	9
Max 300kg Min 2kg e = 0.1kg B.M.L 19.0	TARE CLEAR	0	ENTER

(Wireless functionality optional)


Key Function

- 1. ON/OFF: Power on or power off.
- 2. ZERO: Reset display to 0.0 kg display. Press and hold for 3 seconds to enter device settings.
- 3. M1-5: Saving pre-tare values (up to 5)
- 4. PRE-TARE: Pre-tare the known weight of an object (ex: chair) before beginning measurement.
- 5. TARE: Allows user to deduct weight from reading after measurement
- 6. PRINT: When printer or PC is connected to the scale, press this key to print results
- 7. BMI: Calculation of Body Mass Index
- 8. HOLD: Determine stable weighing value used when weight is unstable. Press and hold for 3 seconds to enter time setting.
- 9. 0-9: For entering digits.



Body Mass Index is a height to weight ratio, and is calculated by the following formula

Definitions Stable symbol: Indicate that weight is stable. Zero symbol: Weight is at zero Negative weight: Weight under zero. Low battery: Battery needs to be charged or replaced.



B.M.J. **6 6.8**

IV. Using Device

A. Basic Operation

Switch on the device using will key. The device will automatically perform self-calibration, displaying software version.

Once "0.00 kg" appears on indicator, device is ready for measurement.

Note: If "0.00 kg" does not display on indicator, press the device.

Guide subject to stand upon the measurement platform. After the weight has stabilized, the "stable" symbol will appear on indicator.

Note: If subject's weight exceeds scale capacity, indicator will display "Err" prompt due to overload.

B. Hold

The hold function determines average weight, designed to be used if subject's weight will not stabilize (ex: an active child).

Note: if fluctuation is too severe, average weight determination will be difficult and hold may not function correctly

- 1. Switch on the device normally.
- IOI D key. "HOLD" will be displayed on the indicator. 2. Press the
- 3. Guide subject to stand on measurement platform.

4. After a few seconds, the average weight will be displayed on the indicator. This weight will be locked - at this point, subject can step off from device.

5. To release the locked weight, press the weight key again to return to the device to normal mode.

Note: Hold function can be activated before or after subject stands on measurement platform. However, if subject finds it difficult to stand still, we recommend activating Hold after subject stands on platform.

C. BMI

1. In normal mode, press the

BM11 kev to enter BMI mode.

2. Display will show last recorded height. Left-most digit will flash.



key to zero

3. Enter height using numeral keys (ex: 170 cm). Input will automatically

move to next digit. Press manually move to next digit.

to confirm.

key to re-input. Press

TARE

kev to

4. After inputting height, press 5. Proceed to weigh subject as usual. Indicator will display weight, height, and BML

NOTE: Hold function can be used at this time if weight is unstable

key to return to normal mode. 6. Press

BMI (w/HM200D or HM201D)

1. Ensure HM200D/HM201D is plugged into indicator.

BM11 2. In normal mode, press the key to enter BMI mode.

3. Proceed to weigh subject as usual. Indicator will display weight, height, and BML

4. Lower stopper on HM200D/HM201D until it touches top of subject's head. Device will automatically calculate BMI based on change in height and weight.

NOTE: Hold function can be used at this time if weight is unstable



key to return to normal mode.

Category	BMI (kg/m ²)	Risk of obesity-related disease
Under	< 18.5	Low
Normal	18.5-24.9	Average
Over	24.9-29.9	Slightly Increased
Obese I	30.0-34.9	Increased
Obese II	35.0-39.9	High
Obese III	> 40	Very High

(World Health Organization adult BMI standards)

D. Tare

5. Press

The tare function allows the user to deduct the weight of objects from the device's measurement result.

1. Place object that needs to be tared onto measurement platform.

TARE

2. Press key after stable symbol appears on indicator. Display will indicate "0.00 kg".

3. Guide subject (plus tared object) to be weighed upon measurement platform. Conduct measurement.

4. To clear tare value, remove all objects from measurement platform,

and press

tare key.

E. Pre-Tare

The Pre-Tare function is used to subtract the known weight of a substance prior to weighing. The device can store 5 sets of pre-tare values.

Pre-tare values can be stored using two different methods: "Load Weight", or "Input Manually".

After pre-tare weights have been stored, they can be recalled by holding

the key for 3 seconds.

С лике С лике иза истират С лике иза истират С лике иза истират Вт рямт ам ноло 2 лике и 2 лике и - 3 4 5 6 и - 3 4 0

A. Load Weight




Press Press key to store pre-tare weight; the indicator will make a beep sound.

0-0-	rîn 🔊	877	PRINT	BMI	HOLD
	ב וו ג ע און ב	ZERO	1	2	3
ſ		MI-5	4	-5	6
	cm	PRE- TARE	7	8	9
		TARE	CLEAR	0	ENTER


B. Input Manually

· · · · ·	
DESCRIPTION	EXAMPLE
Press key. Left-most digit will begin blinking. If no further action is taken within 6 seconds, indicator will return to normal mode	C NVM THE HAD WE FROM THE PRINT BY HOLD C THE
While digit is blinking:	
Enter pre-tare weight using 0~9 keys.	
Ex: to pre-tare 5.0 kg of weight, press 0-0-5-0.	TARE CLEAR O ENTER
Ex: to pre-tare 13.5 kg of weight, press 0-1-3-5.	
Press key to confirm the pre-tare weight.	
Indicator will display minus sign to the left of pre-tare weight value.	



C. Recall Pre-Tare Weight

DESCRIPTION
Press and hold key for 3 seconds. Indicator will display pre-tare value M1 first. The pre-tare value will flash.



NOTE: Pre-tare weight must be under max capacity, otherwise screen will

show 0.00 after key is pressed, and the operator will have to re-input pre-tare settings.

F. Print

If thermal printer is connected to indicator, results can be printed by

pressing key.

USB Port USB Port Port Power socket

NOTE: Thermal printer needs to be powered by adapter

V. Device Setup

A. Setting Time & Date



Press and hold key for 3 seconds to enter Time Setting mode.

Example: Inputting 2008, Dec 25, 8:00am

	Year Setting
מחחב	Enter year using numeral keys 0-9.
2008	(The second sec
	Press key once completed to
	proceed to month & date setting.
	Month & Day Setting
	Enter month, followed by day using
12,25	numeral keys 0-9.
	Ex: December 25th is "12.25".
	Input 1-2-2-5.
	HOLD
	Press E
	proceed to time setting.
	Time Setting
	Enter time (24hr format) using
	numeral keys 0-9.
	Ex: 08:00am is input by pressing
	0-8-0-0.
	HOLD
	Press key once completed to
	confirm time settings and proceed to confirmation.
	Device will display new time and
	date settings, cycling between year,
2008 1225 0800	month & day, and time.
	, , , , , , , , , , , , , , , , , , ,
	YYYY→MM.DD→:HH:MM
	Proce HOLD Kow to return to
	Press key to return to normal weighing mode.
	normal weighing mode.

B. Device Setup

When the device is switched on, press and hold the 2 key for about 3 seconds, until the display shows the "SETUP", followed by "A.OFF" (first option in setting menu).

In device setup menu:

TARE 2510 IOLD

to togale next menu option

to toggle previous menu option

to confirm selection / enter submenu

Auto Power-Off: Instruct device to shut off automatically after a certain period of time.

Press

to togale between options (120 sec / 180 sec / 240 sec / 300 TARE to confirm selection.

sec / off), and



Buzzer/Beep:

When function is turned on, beeping noise will be made when: indicator is turned on, keys are pressed, and weight is stable.



to toggle between on/off, and

TARE

kev to confirm

selection.



Hold Stop: When Hold Stop is "on", Hold will deactivate after subject leaves measurement platform.

Press

to toggle between on/off, and



key to confirm selection.

Language: Set thermal printer language

HOLD to toggle between English, Italian and Polish. Press Press L key to confirm selection.





Font size: Set thermal printer font size.

to toggle between normal and double (larger). Press Press key to confirm selection.





Bluetooth (optional): If device has Bluetooth module installed, Bluetooth function can be turned on or off.

Press to toggle between on/off, and

TARE

to confirm selection.

Wi-Fi (optional): If device has Wi-Fi module installed, Wi-Fi function can be turned on or off.

Press



to toggle between on/off, and **to** confirm selection.

BPSEF

Wi-Fi Setting (optional): If device has Wi-Fi module installed, this option will appear.

to toggle between "Auto" and "PKEY". Press Press confirm selection.

If "Auto" is selected, weight measurement will be automatically sent to connected printer or device. If "PKEY" is selected, transfer will occur

manually only after

key is pressed.

31

VI. Setup USB Connection to PC

For successful connection, PC hardware connected to device must be compatible with USB 2.0 or above. Operators should select a USB cable length that is most suitable to the operating environment.

1. Charder Smart Data Manager can be used to connect the device to a PC. The software program can be downloaded from the Charder website:

[LINK URL] https://www.chardermedical.com/download.htm

2. Connect USB cable to device indicator and PC. Follow installation instructions.

Program Setup

1. After installation of Charder Smart Data Manager is complete, software will automatically search for COM port. Press [**Connect**]. Once connected, **[Connect]** button will change to **[Disconnect**].

Gross Weight	0.0	kg	First Name	Enter	
Tare Weight	0.0	kg	Last Name	Enter	
Net Weight	0.0	kg	Patient ID	Enter	
Height	0.0	cm	Date of Birth	31 / 12 / :	L990 🗐
BMI	0.0		Gender	Male F	emale
Data Auto	Mar	nual			
Please press "Conr Update Time: Model:	nect".		Collect	Clear Sa	ave as

Conducting Measurement

1. Input subject's first name, last name, patient ID, date of birth (DD/MM/YYYY), gender, and height (for BMI calculation) into software if needed. Press **[Clear]** to clear all input.

		kg	First Name	Jane	
Tare Weight	0.0	kg	Last Name	Doe	
Net Weight	0.0	kg	Patient ID	20190201	
Height	167.0	cm	Date of Birth	31 / 12	/ 1965 📊
BMI	0.0		Gender	Male	Female
Data A	Nuto Mai	nual			

NOTE: information can also be input after weight measurement.

2. Conduct measurement. If **[Auto]** is selected, results will be transmitted from device to software automatically and displayed on the left of screen. If **[Manual]** is selected, user must press "Collect".

		kg	First Name	Jane
Tare Weight	0.0	kg	Last Name	Doe
Net Weight	72.5	kg	Patient ID	20190201
Height	167.0	cm	Date of Birth	31 / 12 / 1965 🗐
BMI	26.0		Gender	Male Female

Saving & Printing Results

1. Press **[Save as]** to save measurement results as .csv file on PC. Default file name is same as user ID. (ex: 20190201.csv) To track changes and multiple measurements for the same subject, we recommend not changing the default file name.

Gross Weight	72.5	kg	First Name	Jane
Tare Weight	0.0	kg	Last Name	Doe
Net Weight	72.5	kg	Patient ID	20190201
Height	167.0	cm	Date of Birth	31 / 12 / 1965 🗐
BMI	26.0		Gender	Male Female
Data	Auto Ma	nual		
Data update				

2. Result example:

	A	В	С	D	E	F	G	H	Ī	j	
1	Patient ID	First Name	Last Name	Date of Bi	Gender	Gross We	eis Tare We	igł Net Weig	ht Height	BMI	
2	20190201	Jane	Doe	31/12/1965	Male	72.4 kg	0.0 kg	72.4 kg	167.0 cm		26
3											
4											
5											

If previous results were saved in "20190201.csv", new results also need to be saved as "20190201.csv" (overwriting old file) in order to save multiple results for the same subject.

	А	В	С	D	E	F	G	Н	Ι	J
1	Patient ID	First Name	Last Name	Date of Bi	Gender	Gross Weig	Tare Weigł	Net Weight	Height	BMI
2	20190201	Jane	Doe	31/12/1965	Male	72.4 kg	0.0 kg	72.4 kg	167.0 cm	26
3	20190201	Jane	Doe	31/12/1965	Male	75.2 kg	0.0 kg	75.2 kg	167.0 cm	27
4										

Results will be saved in chronological order of measurement.

3. Press the printer icon to print out result using a printer connected to the PC.

調査の						First Name	kg	72.5	Gross Weight
		5			Doe	Last Name	kg	0.0	Tare Weight
0190201 Jane	201	÷	Patient ID First Name	-	201902	Patient ID	kg	72.5	Net Weight
Doe 965	31/12/196 Male	n :	Last Name Date of Birth Gender	/ 1965 🗐	31 / 3	Date of Birth	cm	167.0	Height
.2 kg .0 kg	75.2	-	Gross Weight Tare Weight	Female	Male	Gender		26.0	BMI
.0 cm	167.0 27.0	:	Height BMI				nual	Auto Ma	Data
.2	75.2		Net Weight Height	Temale	Male	Gender	nual		

VII. Wireless Connection

If the device has the wireless or bluetooth module installed, the indicator can transmit measurement results wirelessly. Please see Charder wireless or bluetooth software instructions for details.

VIII. Troubleshooting

Before contacting your local Charder distributor for repair service, we recommend considering the following troubleshooting procedures:

Self-inspection

1. Device will not power on

- If battery power is depleted, replace with new batteries
- If batteries are not used, check if the power adapter is plugged into the device properly. Check if power adapter is plugged into mains properly

2. Indicator showing "0000" ZERO SPAN out of range

- Interference due to factors such as RF disturbance or ground vibration. Relocate device to location without interference and try again
- Unstable platform feet adjust platform feet according to bubble level indication (clockwise to retract, counter-clockwise to extend) and try again
- External objects interfering with measurement platform. Clear platform of objects and try again
- Device may not function properly on soft surfaces such as carpets or lawns. Relocate device to location with solid, stable floor
- If the steps above cannot resolve the problem, re-calibration may be required to correct weighing accuracy

3. Connection failure for data transmission to PC or printer

- Ensure wires are connected correctly between indicator and PC or printer
- Ensure printer is supplied with power. Ensure PC software is set up properly as indicated in this manual

Distributor support required

If the following errors occur, we recommend contacting your local Charder distributor for repair or replacement services:

1. Device will not power on

- Faulty on/off key
- Broken or damaged wires causing short circuit or faulty connection
- Safety fuse burnout
- Faulty adapter

2. Indicator damage

- Possible hardware defects include: uneven brightness in LCD screen, blurred text, smeared rainbow screen, incorrect decimal display
- Unable to save or read data
- Indicator shows "ERRL" after device is switched on
- Keys not responding
- Buzzer malfunction

Error Messages

Error Message	Reason	Action
Lo	Low battery warning Voltage of battery is too low to operate device	Replace batteries, or plug in adapter
Err	Overload Total load exceeds device's maximum capacity	Reduce weight on measurement platform and try again
ErrH	Counting Error (too high) Signal from loadcells too high	Error normally caused by faulty loadcell or wiring. Please contact distributor
ErrL	Counting Error (too low) Signal from loadcells too low	Error normally caused by faulty loadcell or wiring. Please contact distributor
00000	Zero count over calibration zero range +10% while power on	Re-calibration required. Please contact distributor
00000	Zero count under calibration zero range -10% while power on	Re-calibration required. Please contact distributor
Err.P	Program Error Fault with device software	Error normally caused by faulty loadcell or wiring. Please contact distributor

IX. Product Specifications			
Model		MS3910	
Dis	play	DP3710	
	Capacity	300 kg x 0.1 kg	
Weight	Accuracy	± 1.5e	
Measurement	OIML	Class III	
	LCD Screen	1.0-inch LCD screen (5 1/2 digits)	
	Device	340(W) x 540(D) x 970(H) mm	
Dimensions	Platform	340(W) x 450(D) x 90(H) mm	
	Column	850 mm	
Device	Weight	10.2 kg	
Key Functions		On/Off, Zero, Print, BMI, Hold, Pre-Tare, Tare, Clear, Enter, 0~9, M1-5	
Data Transmission		USB, Wireless module (optional) NOTE : Device should be connected to network by qualified distributors only	
Power Supply		Rechargeable battery pack (optional) or 6 AA batteries / Power adapter	
Operation Temperature & Humidity		0 °C ∼40 °C 15% / 85% RH	
Standard Accessories		(see accessory list)	
Optional Accessories		Thermal Printer, Height Meter	

The device is only compatible with the power adapters specified below.

AMP VOLTAGE	DRAWING NO.	CE APPROVED TYPE NO. / MODEL NO.	ТҮРЕ	Adapter plug
12V 2A	AD-8058(AD-0521)	UE24WU-120200SPA	US	
	AD-8057(AD-0520)	UE24WV-120200SPA	EU	90 - degree
	AD-8056(AD-0519)	UE24WB-120200SPA	UK	
	AD-8074(AD-0534)	UE24W4-120200SPAS	AU	

Standard Accessories

No.	Accessories	Item	Spec.	Qty.
1		USB cable	WR-4001	1
2		12V Adapter		1
3	He bis (phone)	User manual	CD-IN-00327	1
4	E.	Screw (Castor wheel version)	M4*20	4

X. Declaration of Conformity

This product has been manufactured in accordance with the harmonized European standards, following the provisions of the below stated directives:

CE 2460	93/42/EEC as amended by 2007/47/EC Medical Device Directive	
CE M year	2014/31/EU Non-automatic Weighing Instruments Directive	

Please see separate document showing on sticker of device for above CE marking.

Authorized EU Representative:



Obelis s.a.

Bd Général Wahis, 53 B-1030 Brussels Belgium



Manufactured by: Charder Electronic Co., Ltd. No.103, Guozhong Rd., Dali Dist., Taichung City, 41262 Taiwan (R.O.C.)

CD-IN-00327 REV 005 04/2021